Cell Viability/Recovery Specimen Shipment Log

Purpose: To record information about contents of cell viability specimen shipment and receipt of specimens at the JHU Cell Viability Laboratory located at the Johns Hopkins University.

When: Ship on Monday, Tuesday, or Wednesday. Avoid shipments 2 days prior to weekends and holidays. (Each clinic is assigned a month to send in a QC specimen for cell viability tests.)

By whom: Clinic coordinator or laboratory personnel responsible for freezing and shipping of specimens.

Instructions to shipper:

- Process whole blood for leukocytes for (aliquot #7) for cell viability purposes according to LSOCA cell separation protocols.
- Freeze aliquot #7 in a Mr. Frosty or comparable slow freeze at 1° C per minute in a 70° C freezer for a minimum of 4 hours.
- Immediately ship aliquot #7 on the same day specimen was processed with dry ice. Ship on Monday, Tuesday, and Wednesday only. Use the Federal Express air bill provided by Thermo Fisher a third party billing number is preprinted on the air bill.
- Ship via Next Day Arrival using the specimen shippers provided by Thermo Fisher to:

Stacey Meyerer

Johns Hopkins University Bloomberg School of Public Health 615 North Wolfe Street, Room W6614 Baltimore, MD 21205 (410) 955-9810

On the day of shipment:

- Fax the completed Cell Viability/Recovery Specimen Shipment log (VC) to the JHU Lab at (410) 614-2640 to serve as a SHIPPING ALERT.
- Send a shipping alert e-mail to Stacey Meyerer (smeyerer@jhsph.edu), Shivam Chandan (schandan@jhsph.edu) and Akila Hadji (ahadji@jhsph.edu) with the tracking number and clinic/lab contact information.
- Include a hard copy of this specimen log in the box with the specimens; keep a copy of this form for your files.

1. Clinic ID code:				4. Name of person preparing shipment:		
				rvaine of person prepar	ing simplificate.	
2. Date specimen(s) shipped:						
. 11	_	day month y	/ear			
		e specimens sent				
(# of leukoc	yte aliquots)	:				
ecimen shipme	ent informa	tion				
_			! in items 5 thru 7.	. Indicate number of sp	ecimen tube(s) containii	
eukocytes.	J	1 11		J 1		
a.	b.	c.	d.	e.	f.	
		Specimen	# of			
Patient ID#	Visit ID	collection date (day-month-year)	leukocyte aliquots	Date shipped	Initials	
ratient iD#	VISIL ID	(day-month-year)	anquots	Date shipped	muais	
					-	
nments (to be	completed l	oy clinical center)				
innents (to be	completed i					